

BREWSTER ANGLE

Aim: To find the refractive index of transparent material by measuring Brewster angle.

Apparatus: Breadboard, Laser diode, polariser rotator, plate of transparent material (microscope glass slide), rotation stage, photodetector, Detector output unit.


Theory: A beam of light incident on a dielectric transparent material can be resolved into parallel (P) and orthogonal (S) components. These components have different reflection coefficient and Brewster discovered that at a particular angle of incidence θ_B (called Brewster angle), reflection coefficient of P-component goes to zero. At this angle direction of reflected and transmitted beam are orthogonal to each other.

By Snell's law,

$$\tan \theta_B = n,$$

n is the refractive index of the material.

Procedure:



- a) Mount diode laser to the laser mount.
- b) Switch on the laser and place the polariser rotator & analyser in front of it so as to make the E field parallel to the breadboard.
- c) Mount the glass slide on rotation stage.
- d) Orient the microscope slide to reflect the laser beam back into the laser output aperture.
- e) Rotate the glass slide slowly and note the corresponding degree with intensity of the reflected beam from the glass slide.
- f) The intensity has a minimum (almost zero) at Brewster angle θ_B .
Knowing θ_B refractive index is calculated from the formulae, $n = \tan \theta_B$

• **Result:** Calculated the refractive index of transparent material (glass slide) by measuring Brewster angle.

Note:

Initially mount the glass slide, note that the Rotational stage is at zero Degree & the glass slide is parallel to the optical axis. Switch on laser & place a Screen in front of Glass slide so as we can observe the transmitted beam.

Now rotate the Rotational stage with glass slide at keep it at around 53° - 56° and observe transmitted & reflected beam.

Carefully rotate the polariser analyser & rotator by observing the two beams. It can be seen that at a particular degree (rotated polariser analyser & rotator) the reflected beam comes zero whereas transmitted light can still be seen. This is the Brewster angle for that glass plate. This method is used to partially polarise a laser beam.

Only when the beam is partially polarised, you can observe the Brewster angle of a particular material.